

REF: AGX-5-000-0882

LOT: -

Product Identity: C-terminus of Heat-labile enterotoxin B chain (aa 168-319)

Clostridium perfringens

Source: *E.coli*

Size: 19.2 kDa

Description: Protein is expressed with His-tag (N-terminally) and purified with affinity chromatographic techniques

Purity: $\geq 98\%$ by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: 20 mM Tris (pH 8.0), 300 mM NaCl

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only

Avoid repeated freeze-thaw cycles