

REF: AGX-5-000-0713

LOT: -

Product Identity: 15 kDa lipoprotein (aa 12-142)
Treponema pallidum (strain Nichols)

Source: *E.coli*

Size: 16.8 kDa

Description: Protein is expressed with His-tag (N-terminally) and purified with affinity chromatographic techniques

Purity: $\geq 98\%$ by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: 20 mM Tris (pH 8.0), 300 mM NaCl

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only
Avoid repeated freeze-thaw cycles