

**REF:** AGX-5-000-0824

**LOT:** -

**Product Identity:** Full length VCA-p18 (aa 1-176), (ORF: BFRF3)  
Epstein-Barr-Virus (strain B95-8) (HHV-4)

**Source:** *E.coli*

**Size:** 18.9 kDa

**Description:** Protein is expressed with His-tag (C-terminally) and purified with affinity chromatographic techniques

**Purity:**  $\geq 70\%$  by 12.5% SDS-PAGE

**Appearance:** clear solution

**Quantification:** BCA Protein Assay

**Concentration:** Lot specific

**Buffer:** 8 M Urea, 100 mM NaH<sub>2</sub>PO<sub>4</sub>, 10 mM Tris

**Shelf Life:** 3 years at -70°C or below

**Note:** This product is intended for research use only  
Avoid repeated freeze-thaw cycles