

**Article:** AGX-5-000-0947

**LOT:** -

**Product Identity:** Early antigen EA-p54 (aa 306-390 dual hybrid)  
Epstein-Barr-Virus (strain B95-8) (HHV-4)

**Source:** *E.coli*

**Size:** 28 kDa

**Description:** Protein is expressed with His-tag (N-terminally) and purified with affinity chromatographic techniques

**Purity:**  $\geq 80\%$  by 12.5% SDS-PAGE

**Appearance:** clear solution

**Quantification:** BCA Protein Assay

**Concentration:** Lot specific

**Buffer:** 20 mM Tris (pH 8.0), 300 mM NaCl

**Shelf Life:** 3 years at -70°C or below

**Note:** This product is intended for research use only  
Avoid repeated freeze-thaw cycles