

Article: AGX-5-000-0959

LOT: -

Product Identity: EBNA-1 (aa 408-641)
Epstein-Barr-Virus (strain B95-8) (HHV-4)

Source: *E.coli*

Size: 36 kDa

Description: Protein is expressed with His-tag (N-terminally) and purified with affinity chromatographic techniques

Purity: $\geq 70\%$ by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: 20 mM Tris (pH 8.0), 300 mM NaCl

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only
Avoid repeated freeze-thaw cycles