

**REF:**     AGX-5-000-0703

**LOT:**     -

**Product Identity:**    Unique N-terminus of Capsid protein VP1 (aa 1-227)  
                              Human Parvovirus B19

**Source:**                *E.coli*

**Size:**                    26.2 kDa

**Description:**         Protein is expressed with His-tag (C-terminally) and purified  
                              with affinity chromatographic techniques

**Purity:**                 ≥ 98% by 12.5% SDS-PAGE

**Appearance:**         clear solution

**Quantification:**      BCA Protein Assay

**Concentration:**        Lot specific

**Buffer:**                20 mM Tris (pH 8.0), 300 mM NaCl

**Shelf Life:**            3 years at -70°C or below

**Note:**                  This product is intended for research use only  
                              Avoid repeated freeze-thaw cycles