

REF: AGX-5-000-0864

LOT: -

Product Identity: C-terminus of Non-capsid protein NS-1 (aa 403-671)
Human Parvovirus B19

Source: *E.coli*

Size: 31.5 kDa

Description: Protein is expressed with His-tag (N-terminally) and purified with affinity chromatographic techniques

Purity: $\geq 70\%$ by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: 8 M Urea, 100 mM NaH₂PO₄, 10 mM Tris

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only
Avoid repeated freeze-thaw cycles