

Quotation / Order Antibody Service

Seramun Diagnostica GmbH
 Spreenhagener Str. 1
 15754 Heidensee, Germany

T +49 33767 791-10
orders@seramun.com

| Quotation Request | Order |
|------------------------------|--------------------------------------|
| Customer | Contact Person |
| Company / Institute | Name |
| Department | Telephone |
| Street, Number | Email |
| Postcode, City | Quotation no. (if available) |

| Immunization order | Number of animals per antigen | Additional services* |
|---------------------------------------|-------------------------------|---|
| Polyclonal antibodies Sheep | | IgG preparation from antiserum |
| Rabbit | | IgG preparation from antiserum |
| Chicken | | IgY preparation from egg yolk |
| Monoclonal antibodies Mouse | | IgG preparation from cell culture supernatant |

*We offer a broad range of additional services (e.g. antibody conjugation). Please contact us for an individual quotation.

| | | | |
|--|-----------|------------------------------------|--------|
| Antigen Information | | | |
| Name: | | | |
| Size: (MW in kDa without tag or number of amino acids for synthetic peptides) | | | |
| Origin: | | | |
| Native | | | |
| Recombinant | His-tag | GST-tag | Other: |
| Synthetic Peptide | Sequence: | | |
| Concentration: | Volume: | Buffer, additives (if applicable): | |
| Storage temperature: | | | |

Antigen requirements:

- Non-toxic (toxin antigens must be inactivated)
- Non-infectious (infectious antigens must be inactivated)
- Presence in physiological buffer (e.g. PBS, TBS), ± pH neutral (6-8)
- NaCl ≤ 0,3 M
- Imidazole ≤ 0,2 M
- Detergents < 0,2 %
- Urea ≤ 1,0 M
- EDTA < 0,001 M
- TRIS ≤ 0,1 M
- No preservatives (e.g. sodium azide)

Date, signature: